

MedVision



Specification

Power Source	AC 220V, 50/60Hz
Power Consumption	210W
Image Sensor	Sony 1/3" Chip
Definition	410,000 pixels
Min Illumination	0.1 Lux
Video Output	Composite : 1Vp-p
Halogen Lamp	150W(15V)
Light Intensity	Adjustable(0-150W)
Memory Board	1, 2, 4 Capture (NTSC/PAL)
Dimension	300(W) x 285(D) x 125(H)
Weight	Approx. 10 Kg

Standard Set

Monitor (LCD or CRT)	1 EA
Camera Gun	1 EA
Cart	1 EA
Otoscope(0°, Ø4, 50mm)	1 EA

Cart Size

CRT Cart	439 x 450 x 1286.6
LCD Cart	529.4 x 450 x 1254

- Remarkably Much Higher Quality, Very Competitive Price, Long Life Sustainable Operation
- Compact & Smart Design, High Definitive Image Display, Easy Grip of Gun Type Cable

MEDSTAR's Attaining CERTIFICATES



MEDSTAR